

Swissmeda Cloud Imaging requirements

Highlights

- The workstation or tablet must be connected to Internet
- The minimum internet connection to review 3D X-ray is 50 Mbps down link / 20 Mbps up link
- Each User must have a Carestream Dental account created via an invitation by your Admin user
- Each User must be logged with the Carestream Dental account
- Latest browser version of Microsoft Edge or Google Chrome are recommended for the Web viewer
- Optimal 3D rendering requires minimum CPU & RAM and specific browser configuration

See below more details that you may want to review with your IT professional.

Swissmeda Cloud Imaging Viewer Workstation requirements

• Latest versions of Google chrome, MS Edge, Firefox, Apple Safari

Note: specific configuration may be required see Configure / Troubleshoot rendering problems section

- 3D Image viewing requirements
 - o Intel i5 processor or higher / AMD Ryzen 7 With Radeon or higher
 - 8 GB RAM minimum
 - o Ideally Graphics card

Workstation Internet Connection

Recommended for less than one minute download time of one exam

- 20 Mbps down / 10 Mbps up for 2D X-ray images
- 50 Mbps down / 20 Mbps up for CBCT (e.g: 5x5 to 8X9 standard resolution)
- 100Mbps down / 50 Mbps up for CBCT (e.g: 8x9 to 12x10 standard resolution)
- 200Mbps+ down / 100 Mbps up for CBCT (e.g: 12x10 to 16x17 standard resolution)

<u>Click here</u> to check the internet from each of your site from an acquisition workstation, from the CS Imaging server and from a review workstation.

Swissmeda Cloud Imaging Viewer Tablets requirements

• iPad Pro, Galaxy Tab S8+ with the latest Android or iOS version installed. Note: Only 2D images and small size 3D images (up to 20 Mb) can be reviewed

User requirements

- Customer must designate a portal administrator who invites and manages user accounts
- Each user must be always logged to Swissmeda Cloud Imaging Viewer or to CS Imaging

CS Imaging Workstation requirements

- CSI Suite 3.11.0.24 / CS Imaging 8.0.24 or above
- Actively supported version of Windows 11 or Windows 10
- See CS Imaging System Requirements at <u>click here</u>

A Service contract for Carestream Dental equipment, imaging software and DPMS



Configure / Troubleshoot rendering problems

- Make sure your graphics driver and browser are up to date
- In case of unresponsive Page or lack of performance of the viewer, see below
- if you have a Graphic Card (GPU), you may need to enable high performance for better rendering, see below

In case of Unresponsive Page or lack of performance of the Viewer in Chrome or Edge

Open Google Chrome or Microsoft Edge and in the address bar, enter:

- chrome://flags/#use-angle or edge://flags/#use-angle
- select D3D11onD12
- Relaunch button at bottom-right to apply the changes.

In case of Unresponsive Page or lack of performance of the Viewer in Firefox

- Open Firefox and in the address bar, enter "about:config".
- In the preference's search bar, enter "webgl.angle.forced3d11".
- Set the option to "true".

Only if you have a workstation with a Graphic card (GPU) for a better rendering

Open Google Chrome or Microsoft Edge and in the address bar, enter:

- chrome://flags/#ignore-gpu-blocklist or Edge://flags/#ignore-gpu-blocklist
- Enable

You may need to Select High-performance GPU for Microsoft Edge or Chrome in the Windows Graphics settings

Windows-Start-icon /Settings/System/Display/Graphics/Edge or Chrome/Options

- Change Graphic preference to High performance for Chrome or Edge
- **Restart** the browser to apply the changes.

Syst	em	>	Disp	olay	3	Grap	ohi	cs	
Defau	lt sett	ing	5						
Change	detault	Iraph	ics setti	ngs					
Custo	m opi	tion	s for	apps					
Add an a	pp								
Deskto	op app					~			
Brow									
Find an a	app in th for it. Yo	ie list u mig	and sel	ect it, th to rest	hen c tart ti	hoose cust he app for	om g your i	raphics change	to
Find an a settings take effe Search	app in th for it. Yo ct. this list	ne list nu miç	and sel	ect it, th i to rest	hen c tart ti	hoose cust he app for	om g your i	raphics	s to
Find an a settings take effe Search Filter by	app in th for it. Yo ct. this list Apps I	e list u mic	and sel ht need	ect it, th s to rest	hen c tart ti	hoose cust he app for	om g your i	raphics change	s to
Find an a settings take effer	app in th for it. Yo ct. Apps I Google (High per C. Progr	confi confi confi chron	and sel iht need gured the snce	ect it, th i to rest	hen ci tart ti	Appleration	om g your i	raphics change	s to
Find an a settings take effe	app in th for it. Yo ct Apps I Google (High per CAProgr	confi Chron	and sel ht need gured ne ence	ect it, ti i to rest	hen clart ti	hoose cust te app for C	om g your i	raphics change	s to